





Shortcuts in DTX Studio™ Clinic 1.9

Patient Module



■ Patient Module
■ Image Preview



Shortcuts in DTX Studio™ Clinic 1.9





Open Patient Record

- Restore down maximized window
- Finish Spotlight action
- Close dialog or wizard

Switch workspaces:

- | | |
|------------------------------|-----------------------|
| F2: 3D patient / Dental scan | F7: Ceph |
| F3: 3D panoramic | F8: Clinical pictures |
| F4: Tooth | F10: Intraoral camera |
| F5: Panoramic | F12: Recent images |
| F6: Intraoral | |



- Ctrl + Z: Undo
- Shift + Z: Zoom out

Capture with a connected sensor

Ctrl + S: Save (diagnosis)

Ctrl + Y: Redo

Ctrl + M: Maximize / Restore down

Adjust brightness and contrast

Navigate to previous maximized image

Rotate 3D view

Navigate to next maximized image

Capture with intraoral camera (when the livestream is active)

