

DTX Studio™ Clinic 1.8.12 — DTX Studio™ Core 3.9.11 Supported PMS integration — North America

Last updated 16/02/2021

Region	PMS Company	Integration Method
US	Curve Dental	OPP
	Dentrix (G6)	Custom (2-way)
	Dentrix Smart Image (G7)	OPPWeb
	Dentrix Enterprise	Custom
	Eaglesoft (version 20+)	Custom (2-way)
	IMS	OPP
	Mac Practice	OPP
	Open Dental	OPP
Canada	ClearDent	OPP

Didn't you find the Practice Management Software you want to connect with?
Contact us at support@dtxstudio.com.